

### 5.5.13.1 Support Equipment

This Standard applies to all devices that may not be the actual test instrument, but are necessary to support laboratory operations. These include, but are not limited to: balances, ovens, refrigerators, freezers, incubators, water baths, temperature measuring devices (including thermometers and thermistors), thermal/pressure sample preparation devices and volumetric dispensing devices (such as Eppendorf® or automatic dilutor/dispensing devices).

a) All support equipment shall be maintained in proper working order. The records of all repair and maintenance activities, including service calls, shall be kept.

**Deleted:** if quantitative results are dependent on their accuracy, as in standard preparation and dispensing or dilution into a specified volume

~~b) Raw data records shall be retained to document equipment performance.~~

**Formatted:** Bullets and Numbering

~~c) On each day the equipment is used, balances, ovens, refrigerators, freezers, incubators and water baths shall be checked and documented. The acceptability for use or continued use shall be according to the needs of the analysis or application for which the equipment is being used.~~

**Formatted:** Bullets and Numbering

~~d) Temperature measuring devices shall be calibrated or verified at least annually. Calibration or verification shall be performed using a recognized National Metrology Institute traceable reference, such as NIST, when available.~~

**Formatted:** Bullets and Numbering

~~i) If the temperature measuring device is used over a range of 10°C or less, then a single point verification within the range of use is acceptable;~~

**Formatted:** Numbered + Level: 1 + Numbering Style: i, ii, iii, ... + Start at: 1 + Alignment: Left + Aligned at: 1" + Tab after: 1.5" + Indent at: 1.5"

~~ii) If the temperature measuring device is used over a range of greater than 10°C, then the verification must bracket the range of use.~~

e) All other support equipment shall be calibrated or verified at least annually, using a recognized National Metrology Institute, such as NIST, traceable references when available, bracketing the range of use.

**Formatted:** Bullets and Numbering

~~f) The results of such calibration or verification, whether daily or annually, shall be within the specifications required of the application for which this equipment is used or:~~

**Deleted:**

i) the equipment shall be removed from service until repaired; or

ii) the laboratory shall maintain records of established correction factors to correct all measurements.

~~g) If quantitative results are dependent on their accuracy, such as in standard preparation or dispensing or dilution into a specified volume, volumetric dispensing devices (except Class A glassware and Glass microliter syringes) shall be checked for accuracy on a quarterly basis.~~

**Deleted:** ¶  
<#>Raw data records shall be retained to document equipment performance.¶  
¶  
<#>On each day the equipment is used, balances, ovens, refrigerators, freezers and water baths shall be checked and documented. The acceptability for use or continued use shall be according to the needs of the analysis or application for which the equipment is being used.¶

**Formatted:** Bullets and Numbering

**Deleted:** Volumetric